

Advantages of HIDREX Devices

- Operation via App with language control
- Support options via App for therapy optimization
- Recognition of your optimal treatment parameters
- Fully automatic therapy procedure
- Effective current types
- Variable pulse current setting
- Powerful devices, up to 60V and 30 mA
- High safety through protective circuits
- ASE anti shock electronic
- Protection against over-treatment
- Immersion monitoring
- Operation without assistance
- Timer function
- Intelligent microprocessor control



H DREX [®] MEDICAL AND PEHAB



HDREX [®] MEDICAL AND REHAB independent living

About us

Over 30 years ago, in cooperation with the largest dermatology clinics in Germany, HIDREX developed the first iontophoresis device for the treatment of hyperhidrosis. It was introduced into the market in 1987, and has been successful ever since. The HIDREX-therapy is now a worldwide scientifically recognized treatment concept and is known as tap water iontophoresis (TWI). We use the latest technologies to manufacture certified medical devices for innovative, natural healing processes.

- Licensed to manufacture medical products in accordance with the EU Directive
- Certified according to EN ISO 13485
- FDA registered

Distributor:



HIDREX GmbH Otto-Hahn-Str. 12 | 42579 Heiligenhaus | Germany Fon: +49 (0) 20 56 / 98 11 -0 | Fax: -31 info@hidrex.de





TWI Treatment for Hyperhidrosis




www.hidrex.com +49 (0) 20 56 / 98 11 10



Hyperhidrosis

The term hyperhidrosis refers to excessive sweating that goes beyond the body's need for thermal regulation. People with hyperhidrosis sweat excessively, regardless of physical exertion or ambient temperature.

This dysfunction can lead to significant mental, social, and occupational limitations and can thus cause a marked reduction in the quality of life. Almost 5% of the population suffer from focal hyperhidrosis. This disease is wrongly considered by many to be untreatable.



Tap water iontophoresis (TWI)

Iontophoresis is a treatment method that has been scientifically recognized for many years and can result in normally dry skin without medication or surgery.

In comparison to many other therapeutic methods, iontophoresis has no negative long-term effects. The procedure is easy to perform, safe, inexpensive, and can be used even on children. In addition, beyond the purchase of an iontophoresis device, there are no further follow-up costs, so that the treatment is very affordable over the longer term.

Buy HIDREX TWI Devices

For more information visit

our website www.hidrex.com

HIDREX has distributors all over the world. Find your local supplier on hidrex.com/buy. For informations about rental options, rent-to-own programs or money back guarantee contact your competent supplier.

Many insurance companies reimburse the cost of the device.

BUY DEVICE



HIDREX iontophoresis treatment

In HIDREX-therapy, a mild DC or pulse current is passed through the affected body zones. The flow of this electrical current is ensured by a water bath and leads to a sweat reduction on the areas being treated.

The HIDREX treatment is carried out by hand and foot baths or by treatment of the armpits, face, neck, or back with special sponge applicators.

During the initial treatment phase, you should, if possible, treat daily for 15-minutes. After about 10 treatments, your sweat levels should normalise. Afterwards, 1-3 treatments a week are usually sufficient to maintain the results. The therapeutic effect is reversible, so if treatments stop, the sweating will resume after some time. Therefore, the purchase of a home therapy device is recommended.



H DREX ^{® MEDICAL} AND REHAB