



Therapy hands



Therapy feet



Therapy armpits

**HIDREX**® MEDICAL AND REHAB  
independent living



### About us

Over 30 years ago, in cooperation with the largest dermatology clinics in Germany, HIDREX developed the first iontophoresis device for the treatment of hyperhidrosis. It was introduced into the market in 1987, and has been successful ever since. The HIDREX-therapy is now a worldwide scientifically recognized treatment concept and is known as tap water iontophoresis (TWI).

### Advantages of HIDREX Devices

- Operation via App
- 98% customer satisfaction
- Non-surgical and drug free treatment
- Direct or pulsed current for user customization and comfort
- Fully automatic therapy procedure
- Powerful devices, up to 60V and 30 mA
- ASE – anti shock electronic
- Protection against over-treatment
- Immersion monitoring
- 100% made in Germany

#### DEVICE FEATURES



Ionto-EN\_Trifold-D\_07-2022



**HIDREX** GmbH

Otto-Hahn-Str. 12 | 42579 Heiligenhaus | Germany  
Fon: +49 (0) 20 56 / 98 11 -0 | Fax: -31  
info@hidrex.de

# STOP SWEATING!

## TWI Treatment for Hyperhidrosis



#### DEVICE INFO



**HIDREX**® MEDICAL AND REHAB



www.hidrex.com  
+49 (0) 20 56 / 98 11 10



**HIDREX**® MEDICAL AND REHAB



### Hyperhidrosis

The term hyperhidrosis refers to excessive sweating that goes beyond the body's need for thermal regulation. People with hyperhidrosis sweat excessively, regardless of physical exertion or ambient temperature.

This dysfunction can lead to significant mental, social, and occupational limitations and can thus cause a marked reduction in the quality of life. About 5% of the population suffer from focal hyperhidrosis, such as moist hands, feet and armpits.

### Tap water iontophoresis (TWI)

Iontophoresis is a treatment method that has been scientifically recognized for many years and can result in normally dry skin without medication or surgery.

In comparison to many other therapeutic methods, iontophoresis has no negative long-term effects. The procedure is easy to perform, safe, inexpensive, and can be used even on teenager and children.

Convince yourself of the effectiveness of tap water iontophoresis!

### Buy HIDREX TWI Devices

HIDREX has distributors all over the world. Find your local supplier on [hidrex.com/buy](http://hidrex.com/buy). For informations about rental options, rent-to-own programs or money back guarantee contact your competent supplier or for Germany HIDREX directly.

**i** Many insurance companies reimburse the cost of the device.

### HIDREX iontophoresis treatment

In HIDREX-therapy, a mild DC or pulse current is passed through the affected body zones. The flow of this electrical current is ensured by a water bath and leads to a sweat reduction on the areas being treated.

The HIDREX treatment is carried out by hand and foot baths or by treatment of the armpits, face, neck, or back with special sponge applicators.

During the initial treatment phase, you should, if possible, treat daily for 15-minutes. After about 10 treatments, your sweat levels should normalise. Afterwards, 1-2 treatments a week are usually sufficient to maintain the results. The therapeutic effect is reversible, so if treatments stop, the sweating will resume after some time. Therefore, the purchase of a home therapy device is recommended.



CE0633



BUY DEVICE



**STOP SWEATING!**

For more information visit our website [www.hidrex.com](http://www.hidrex.com)

**HIDREX**® MEDICAL AND REHAB